Types of Specimen Collection Containers 收集標本容器分類



Other Specimen Collection Containers 其他收集標本容器		
Containers 容器	Containers Name 容器名稱	Suitable for 適用於
	Sterile Urine Container (Universal container)	Urine, semen, sputum, body fluid, fine needle aspirates, cutaneous samples for: • Bacterial culture tests • Fungal culture tests • STD DNA tests
	Non-sterile Urine Container	 Urinalysis (Urine Routine) Kidney stone analysis
	Early Morning Urine (EMU) Container	Special tests
	24-hrs Urine Container *please refer to "Specimen Collection and Handling Guidelines" section for urine preservatives information	• 24-hrs urine tests
1	Sterile Stool Container	Stool cultureFilmarray Gastrointestinal Panel
94	Non-sterile Stool Container	Stool routineStool Occult Blood
Red cap	Swab (with transport medium)	Wound, eye, female urogenital, throat, rectal samples for: • Bacterial culture tests • STD DNA tests
Orange cap	Swab (with transport medium)	Male urogenital samples for: • Bacterial culture tests • STD DNA tests
Blue cap	NP Swab (with transport medium)	Nasal and nasopharyngeal samples for: Bacterial culture tests Influenza A/B
Black cap	Dry Swab	Wound, eye, female urogenital, throat, rectal samples for: • Bacterial Culture tests • STD DNA tests
Daniel Control	UTM Swab (Universal Transport Medium Swab)	Filmarray Respiratory Panel
	13C-UBT Collection Kit	13C Urea Breath Test for Helicobacter pylori
	CytoRich Vial (Thinprep/ SurePath)	Pap smearHPV Genotyping